



Maintenance

Cleaning

Strap, Spreader Bar Cover, and Harness

- Cleaning Frequency: Facility infection control guidelines.
- Cleaning Agent: Ethyl or isopropyl alcohol cleaning agent.

Overview

- Annual preventative maintenance should be performed by a certified Gorbel Medical technician. Please contact Gorbel Medical to schedule this service.

Before Each Use

- Inspect the following components for damage or excessive wear. If any is found, discontinue using and contact Gorbel Medical.
- **Actuator Strap** - Inspect Strap for fraying, tears, rips, cuts or signs of excessive wear
- **Spreader Bar** - With one hand holding the body of the spreader bar, push and pull on the clip assembly (1 and 2) with the other hand to check for excessive “play.” Spacing between the clip assembly and spreader bar should be virtually imperceptible and the clip assembly should rotate smoothly.

Before Each Use

- **Safegait EMBRACE™ Rehabilitation Harness** - Please refer to the **SafeGait EMBRACE** Instructions for Use (IFU) manual for detailed illustrations to reference component parts listed in the inspection items below. All buckle hardware for Leg Cuff Straps; Donning Straps; and Spreader Bar Strap Connectors should be inspected for cracks or defects;
- Inspect all harness components for fraying, tears or signs of wear (e.g., stitching). Check all Chest and AbBand VELCRO® to ensure it provides intended hold when fastened;
- Check D-Rings for any cracks, broken welds, putting or scaring. Check Spreader Bar Straps for wear or fraying.
- Harness should not be used if any signs of damage or excessive wear are evident.

Annual Maintenance

- Contact Gorbel Medical to schedule a preventive maintenance appointment for the following:
- Inspection and cleaning of the ceiling-mounted rail
- Inspection of the SafeGait 360 Actuator
- Replacement of the Actuator Strap